avnavimant titla	VP ITC CDC CAII										Notes
experiment title	VP-ITC CBS-CAII 10/25/2012										
user	John D. Chodera										
titrant:titrand filename	102512a.itc										
titrant:buffer filename	102512b.itc										
EXPERIMENTAL DETAILS											
temperature estimated Ka		25 C 1.00E+06 1/M									
sample cell volume (approximate)		1.00E+06 1/M									
number of injections		10									
injection volume		10 uL									
dilution factor for single injection (d)		0.9928571429									
dilution factor after final injection (d^n)		0.9308241571									
		desired				actual		error		percent error	
cell concentration c-value		10.00 uM 10.00				10.01 12.01		0.09		0.87% 0.26%	
optimal Rm		5.34				4.98		0.03		0.20%	
syringe concentration		718.29 uM				718.23		8.69	uM	1.21%	
LIGAND											
compound name	4-Carboxybenzenesulfonamide (CBS)										
description	4-sulfonylbenzoic acid at 97% purity						-				Color key
vendor product no.	000044709										fill in before experiment
lot no.	#MKRE3323V										automatically computed
purity	MINICI OCCU	97.00%									automatically computed
molecular weight		201.2 g/mol									
solubility		453 mg/L	2,251.49	uM							
Stock solution preparation											
		desired	typical error		typical percent error	actual		error		percent error	
target compound mass	(most balances need min 10 mg)	10 mg	0.1	ma	1.00%	10.01	ma		mg		balance tared to 50 mL Falcon tube + holder near 0.1 mg uncertainty
buffer volume needed for target mass	(for planning buffer usage)	32.14 mL	0.1	9	1.50%	10.01	9	0.1	9	1.0070	bulance and to be the raison table a holder hear our ing anostrainty
buffer volume needed for actual mass	(use this for actual preparation)	32.17 mL	0.64	mL	2.00%	32.17	mL	0.21	mL	0.65%	via P5000: 4x5 mL + 4.130 mL + 2x4.020 mL
purity-corrected stock solution concentration	(be careful not to exceed solubility)	1,500.00 uM	33.54	uM	2.24%	1,500.12	uM	17.90	uM	1.19%	
Syringe solution preparation											
		desired	typical error		typical percent error	actual		error		percent error	
dilution factor from stock to titrant solution		0.479	0.004		0.74%	0.47878		0.00094		0.20%	
stock solution volume		4.309 mL	0.04		1.00%	4.309	mL	0.012			via P5000: 4.309 mL
buffer volume		4.691 <b>mL</b>	0.05		1.00%	4.691		0.012			via P5000: 4.691 mL
total volume for titrant	(need min 350 uL for VP-ITC)	9 mL	0.06		0.71%		mL	0.017	mL	0.19%	
titrant concentration		718.29 uM	17	uM	2.35%	718	uM	9	uM	1.21%	(this error is included in the 'actual error' in thermodynamic parameters below
PROTEIN											
TROTEIN											
protein name	carbonic anhydrase II (CAII)										
source	Sigma-Aldrich C2522-25mg										
lot no.											
molar extinction coefficient		50070 M-1 cm-1	1								
protein purity (or 100% if unknown)		100.00%									
Dialyzed stock solution concentration											
						actual		error			
absorbance measurement for 1 cm path length	(e.g. NanoDrop)					11.768		0.03			used error from spec manual (for 200-350nm)
concentration of stock solution						235.03		0.60		0.25%	
purity-corrected concentration of stock solution						235.03	uM	0.60	uM	0.25%	
Call calution proparation											Add a section for UV-VIS measurement of protein final dilution
Cell solution preparation		desired	typical error			actual		error		percent error	Add a Section for UV-VIS measurement of protein final dilution
dilution factor		0.04255	0.00058			0.04258		0.00036		0.83%	
stock solution volume		0.511 mL	0.005	mL		0.511		0.004		0.78%	via P1000
buffer volume		11.489 mL	0.11	mL		11.489	mL	0.044	mL	0.38%	via P5000: 2x5 + 1.489
total volume for titrate	(need min 2.1 mL for VP-ITC)	12.00 mL	0.12			12.00		0.044	mL	0.37%	
titrate concentration		10.00 uM	0.14	uM		10.01	uM	0.087	uM	0.87%	(this error is absorbed into the n parameter in the ITC fit)
THERMODYNAMIC PARAMETERS											
THERWIOD TNAMIC PARAMETERS											
		reported from fit	error from fit		percent error from fit	actual		error		percent error	
n (stoichiometry, purity, and V0 correction)		0.915	0.003		0.37%	0.915		0.004		0.42%	Origin fit used purity-corrected protein concentration in cell
K (association constant)		1.20E+06 M-1	2.78E+04	M-1	2.32%	1.20E+06	M-1	3.14E+04	M-1	2.61%	· · · · · · · · · · · · · · · · · · ·
Kd (dissociation constant)		8.33E-07 M	1.93E-08		2.32%	8.33E-07		2.18E-08		2.61%	
DeltaH		-11.27 kcal/mol			0.51%		kcal/mol		kcal/mol	1.31%	
TDeltaS DeltaG		-2.97 kcal/mol -8.30 kcal/mol		kcal/mol kcal/mol	1.98% 0.17%		kcal/mol kcal/mol		kcal/mol kcal/mol	5.00% 0.19%	