

Dialysis ~114  $\mu$ M CAII in PBS

1:2 dilution of 228  $\mu$ M CAII stock w/  
PBS Buffer in Maxi size D-tube dialyzer  
against 1L of PBS buffer w/ stirring @  
4°C (coldroom) O/N.

→ Need to check purity from supplier!

2011 Apr 27 9:36 am

Filtred dialysis buffer

~~Pure~~ Diluted 1 mL "polymerized protein" w/ 1 mL  
filtred dialysis buffer to make ~50  $\mu$ M protein stock

Diluted (1:4) 25  $\mu$ L of protein soln into 75  $\mu$ L filtred  
dialysis buffer ↓

$$A_{280} = 0.71549 \text{ AU}$$

∴ ~50  $\mu$ M protein soln is actually ~~50~~  
57.2  $\mu$ M by spec

2.1 $\mu$ L	40 $\mu$ M	<del>10</del> $\mu$ M	X	MX		10 <sup>6</sup> $\mu$ M / M L
		MX	X	10 <sup>3</sup> $\mu$ M	10 <sup>6</sup> $\mu$ M	52.7 $\mu$ M

Added 1.467  $\mu$ L of 57.2  $\mu$ M protein soln to 0.633  $\mu$ L buffer  
to make ~2.1 mL of ~40  $\mu$ M protein soln.

Actual  $A = 0.48059$  for 1:4 dilution  
of  $\sim 40 \mu\text{M}$  stock

$$A = \epsilon bc$$

$$\frac{0.48059}{50,070 \text{ M}^{-1} \text{ cm}^{-1} (\text{cm})} \times \frac{4}{1} \times 10^3 \times 10^3 = 38.4 \mu\text{M} \text{ protein}$$

Friday @ 2PM (meeting - / Bozer)

ITC

$$[\text{COS}] \quad 1782 \mu\text{M}$$

$$[\text{C41}] \quad 40 \mu\text{M}$$

$$\epsilon_{\text{COS}} = 50,070 \text{ M}^{-1} \text{ cm}^{-1}$$

Made 1:2 dilution of O/N dialysis protein  
stock ( $\sim 228 \mu\text{M}$ ).

$$A = 0.71549 \quad \cdot \quad \frac{4}{1} \times 10^3 \times 10^3 = 57.2 \mu\text{M}$$
$$\frac{50,070 \text{ M}^{-1} (\text{cm})}{\text{cm}}$$

$$57.2 \mu\text{M} \cdot \frac{x}{2.1 \text{ mL}} = 40 \mu\text{M} \quad x = 1.467 \text{ mL of } 57.2 \text{ stock}$$
$$\frac{-2.1 \text{ mL}}{0.633 \text{ mL buffer}}$$

Measured A of  $\sim 40 \mu\text{M}$  solution:

$$A = 0.48059$$

$$c = \frac{0.48059}{50,770 \text{ M}^{-1} \text{ cm}^{-1} (\text{cm})} \cdot \frac{4}{1} \cdot 10^3 \cdot 10^3 = 38.4 \mu\text{M protein}$$

$$\text{[CBS]} \text{ mM} \approx 1.783$$

$$\text{MW} = 201.2 \text{ g/mol}$$

$$\text{vol} \approx 10 \text{ mL}$$

$$\text{~~0.010 L~~ } 0.010 \text{ L} \cdot 1.783 \times 10^3 \frac{\text{mol}}{\text{L}} \cdot \frac{201.2 \text{ g}}{\text{mol}} = 0.0036 \text{ g}$$

Actually measured 0.0034 g CBS

$$\frac{3.4 \text{ mg}}{x} = \frac{3.6 \text{ mg}}{10} \quad x = 9.478 \text{ mL PBS buffer}$$

Note CBS was ~~very~~ static  $\rightarrow$  ~~stuck~~ stuck to weighing paper so had to rinse off w/ buffer.

$$\text{CBS, } \epsilon_{272 \text{ nm}} \approx \text{~~1000 - 1500~~ } 1307 \text{ M}^{-1} \text{ cm}^{-1} \text{ (+126)}$$

$\frac{1}{4}$  dilution of  $\sim 1.783 \text{ mM}$  stock

$A = 0.63718$  @  $272 \text{ nm} \rightarrow$  although  $274 \text{ nm}$  looks like max

$$\frac{0.63718}{1307 (1)} \cdot \frac{4}{1} \cdot 10^3 = 1.950 \text{ mM} \quad \text{actual CBS concentration}$$



Run 1:

\* No control injection

$$[CAII]_{mM} = 38.4$$

instead moved plunger

$$[CBS]_{mM} = 1.783$$

down position once  $\times 10 \mu L$

$$\# \text{ injects} = 10$$

$$\text{Vol injects} = 10 \mu L$$

name = 042711a, saved in DATA/Chodera

run file saved as 10 $\mu L$  - 40 $\mu M$

Run 2:

$$A = \frac{0.48678}{50,000 \text{ M}^{-1} \text{ cm}^{-1} (\text{cm})} \cdot \frac{4}{1} \times 10^3 \times 10^3 = 38.9 \mu M$$

$$[CAII] = 38.9 \mu M$$

$$[CBS] = 1.783 \text{ mM}$$

$$\# \text{ injects} = 15$$

$$\text{Vol injects} = 6.7 \mu L$$

$$\text{name} = 042711b$$

Run 3: Heat of dilution

for run # 2, same params

$$\text{name} = 042711c$$

\* Stopped @ pt # 5 b/c of  
baseline shift after pt # 2

Made 1:2 dilution w/ all remaining protein from dialysis

$$57.2_{\mu\text{M}} \cdot \frac{0.900 \text{ mL}}{x} = 40_{\mu\text{M}} + \frac{60_{\mu\text{L}} \cdot 52_{\mu\text{M}}}{x} = 40_{\mu\text{M}}$$

$$\begin{aligned} x &= 1.287 \\ &- 900 \text{ protein stock} \\ &\hline &0.387 \text{ buffer} \end{aligned}$$

$$\begin{aligned} x &= 85.5 \\ &- 60 \\ &\hline &25.5_{\mu\text{L}} \text{ buffer} \end{aligned}$$

Run #4

Repeat heat of dilution w/ Run #2  
params.

Waved = 042711d

Run #5

$$[\text{CAII}] \approx 40_{\mu\text{M}}$$

$$[\text{CBS}] = 1.783_{\text{mM}}$$

$$\# \text{ injects} = 20$$

$$\text{Vol inject} = 5_{\mu\text{L}}$$

$$\text{name} = 042711e$$

Ran out of protein

added 80  $\mu\text{L}$ , took

Abs after w/ leftover

$$c = \frac{0.39238}{50,070(1_{\text{cm}})} \cdot \frac{1}{1}$$

$$\times 10^3 \times 10^5$$

$$[\text{CAII}] = 31.4_{\mu\text{M}}$$

There was a bubble in the cell, happened during washing and dilution of sample. Data from this run isn't any good. Need to repeat.

5/5/2011

Dilution for UV-Vis spec. CA11

$$38.4 \mu\text{M} \cdot \frac{1}{4} = 9.60 \mu\text{M} \text{ for } 1 \text{ cm PL}$$

for nano-drop PL = 0.1 cm (10mm)

should probably do x10 for nano-drop

$$C = 95.98 \mu\text{M}$$

$$A = \epsilon bc = 50,070 \text{ M}^{-1} \text{ cm}^{-1} (0.1 \text{ cm}) (9.60 \times 10^{-5} \text{ M})$$

$$A = 0.48059 \leftarrow \text{Target Abs btw } 0.1 \text{ and } 1$$

[CA11] stock = 228  $\mu\text{M}$

Will: do 1:2 dilution for dialysis

$$228 \div 2 = 114 \mu\text{M} \text{ after dialysis}$$

§ Sample for nano-drop: (tomorrow)

$$114 \mu\text{M} \cdot \frac{1}{x} = 95.98 \mu\text{M} \quad x = 1.188$$

$$\Rightarrow \text{or } 114 \mu\text{M} \cdot \frac{x}{10 \mu\text{L}} = 95.98 \mu\text{M}$$

$$x = 8.419 \mu\text{L of } 114 \mu\text{M stock}$$

$$1.581 \mu\text{L buff.}$$



## 5/12/11 Heat of Dilution control Experiments

1<sup>st</sup> 80  $\mu$ M CAII

From 99.89  $\mu$ M stock

$$99.89 \mu\text{M} \cdot \underline{1.70 \text{ mL}} = 80 \mu\text{M}$$

$$x = 2.12264 \text{ mL} = 2.123 \text{ mL}$$

— 1.70 mL protein stock

0.423 mL buffer

2<sup>nd</sup> 40  $\mu$ M CAII

Used 80  $\mu$ M from previous run (recovered)

First calc. absorbance w/ nanodrop:

$$\sim 80 \mu\text{M} \rightarrow 75.9 \mu\text{M}$$

(1:10 dilution of  $\sim 80 \mu\text{M} \rightarrow 5 \mu\text{M}$ ?)

After run 051211a  $\rightarrow \sim 75 \mu\text{M}$  CAII diluted to  $\sim 70 \mu\text{M}$ .

To make 40  $\mu$ M stock

$$70 \mu\text{M} \times \frac{x}{2.1 \text{ mL}} = 40 \mu\text{M}$$

$x = 1.17 \text{ mL}$

$$x = 1.70 \text{ mL recovered } 80 \mu\text{M CAII}$$

$$2.1 - 1.7 = 0.4 \text{ mL buffer}$$

SKIP TO NEXT PAGE



\* Had to degas protein sample 2x -  
was very bubbly  
Run name 05/2/11b  
Same params as previous

- 3rd Run CA11 (10  $\mu$ M) 20x Syt injections, 446  $\mu$ M  
CBS - Follow-up to experiments on 5/6/11

Used ~80  $\mu$ M CA11 recovered stock again

$$\rightarrow 70 \mu\text{M} \cdot \frac{x}{2.1 \text{ mL}} = 10 \mu\text{M}$$

$$x = \frac{0.3 \text{ mL stock protein}}{1.8 \text{ mL PBS}}$$

Degassed again

Saved as 05/2/11c



5/19/11

Finish 20 x 5  $\mu$ L injections

1<sup>st</sup> 20 x 5  $\mu$ L CBS into PBS buffer  
300s delay  
240s spacing  
Used 446  $\mu$ M CBS from 5/12/11  
Saved as 051911a

2<sup>nd</sup> 20 x 5  $\mu$ L CBS into ~10  $\mu$ M CA11  
Same as above  
refilled syringe w/ ~150  $\mu$ L remaining  
CBS degassed in small cuvette by  
simply ~~ref~~ moving syringe position back  
to "closed"  $\rightarrow$  purged twice, 10  $\mu$ L  
down injection  
saved as 051911b

3<sup>rd</sup> 20 x 5  $\mu$ L buffer into ~10  $\mu$ M CA11  
Same as above

$\rightarrow$  ran out of protein stocks, had to add more  
buffer to cell (see following page)  
Saved as 051911c

### To do list:

- ~~Protein~~ Protein stock from dialysis O/N
  - store
  - check concentration w/ nano-drop
  - filter + store buffer (dialysis #3 now)

- Check  $\sim 10 \mu\text{M}$  CAII stocks by nano-drop
- Make another  $10 \mu\text{M}$  stock from  $80 \mu\text{M}$  CAII recovered ~~to~~ (same as previous) for buffer into protein run

$$\sim 70 \mu\text{M CAII} \cdot \frac{x}{2.1 \text{ mL}} = 10 \mu\text{M}$$

$$x = 0.3 \text{ mL CAII stock}$$

add  $1.8 \text{ mL}$  PBS buffer

→ or just add to previous  $\sim 10 \mu\text{M}$  CAII stock

$$\sim 70 \mu\text{M} \cdot \frac{0.175 \text{ mL}}{x} = 10 \mu\text{M}$$

$$x = 1.225 \text{ mL } \text{total}$$

$$- 0.175 \text{ mL}$$

added to  $\swarrow$   $1.050 \text{ mL}$  buffer

$\sim 2 \text{ mL}$  of  $10 \mu\text{M}$  CAII stock from 5/12/11

Protein for run #3, 051911c.

Had  $\sim 1.5\text{ mL}$  of  $8.07\mu\text{M}$  stock (abs from nano-drop) prior to addition of extra buffer.

To load needed  $\sim 2.1\text{ mL}$   $\rightarrow 2.1 - 1.5 = 0.6\text{ mL}$  buffer added

$$8.07\mu\text{M} \cdot \frac{1.5\text{ mL}}{2.1\text{ mL}} = 5.76\mu\text{M}$$

That's at the lower end;  
however buffer was added  
to "top off" the cell after  
protein solution was loaded,  
so actual concentration is  
somewhere b/w the two.

\* Can back calculate after run is  
finished b/c no CBS present

Abs was 0.337 for sample taken  
after the titration had finished  $\rightarrow$



Started w/ 1.43ml of [CAII] at unknown concentration. Diluted by 100 $\mu$ l of buffer during the course of titration for a final volume of 1.53 ml

$$[\text{CAII}]_f = 6.73 \mu\text{M}$$

$$6.73 \mu\text{M} \times \frac{1.53 \text{ mL}}{1.43 \text{ mL}} = 7.20 \mu\text{M} \leftarrow [\text{CAII}]_i$$

Starting [CAII] was 7.20  $\mu\text{M}$  prior to titration w/ buffer.