

ITC experiment										Notes
temperature		25	C							
estimated K _a		1.00E+08	1/M							
sample cell volume (approximate)		1.4	mL							
number of injections		10								
injection volume		10	uL							
dilution factor for single injection (d)		0.99285714285714								
dilution factor after final injection (d ⁿ)		0.93082415714980								
		desired				actual		error		percent error
cell concentration		10.00	uM			7.99	uM	0.21	uM	2.67%
c-value		10.00				9.99		0.04		0.42%
optimal R _m		5.34				5.34				
syringe concentration		718.29	uM			698.41	uM	7.94	uM	1.14%
LIGAND										
compound name	CBS	Carboxybenzenesulfonamide (CBS)								
description	4-sulfonylbenzoic acid at 97% purity									Color key
vendor	Sigma-Aldrich									fill in before experiment
product no.	D00944798									fill in during experiment
lot no.	#MKBFE3323V									automatically computed
purity		97.00%								
molecular weight		201.2	g/mol							
solubility		453	mg/L	2,251.49	uM					
Stock solution preparation										
		desired		typical error		typical percent error		actual		error
target compound mass	(most balances need min 10 mg)	10	mg	0.1	mg	1.00%		10.1	mg	0.99%
buffer volume needed for target mass	(for planning buffer usage)	24.11	mL							
buffer volume needed for actual mass	(use this for actual preparation)	24.35	mL	0.49	mL	2.00%		25.099	mL	0.27%
purity-corrected stock solution concentration	(be careful not to exceed solubility)	2,000.00	uM	44.72	uM	2.24%		1,940.03	uM	1.03%
Syringe solution preparation										
		desired		typical error		typical percent error		actual		error
dilution factor from stock to titrant solution		0.370		0.003		0.89%		0.36000		0.00176
stock solution volume		1.851	mL	0.02	mL	1.00%		1.8	mL	0.012
buffer volume		3.149	mL	0.03	mL	1.00%		3.2	mL	0.012
total volume for titrant	(need min 2.1 mL)	5	mL	0.04	mL	0.73%		5	mL	0.017
titrant concentration		718.29	uM	18	uM	2.48%		698	uM	8
PROTEIN										
protein name	carbonic anhydrase II (CAII)									
source	Sigma-Aldrich C2522-25mg									
lot no.										
molar extinction coefficient		50070	M-1 cm-1							
protein purity (or 100% if unknown)		80.00%								
Stock solution preparation										
		desired		typical error				actual		error
Dilution for UV-Vis absorbance measurement		0.0625		0.001				0.06250		0.00086
dilution factor for absorbance measurement		10	uL	0.1	uL	1.00%		10	uL	0.1
volume of stock solution		150	uL	1.5	uL	1.07%		150	uL	1.6
volume of buffer solution		160	uL	1.50	uL	1.00%		160	uL	1.603
total volume of dilution										
absorbance measurement								0.66844	A	0.01
concentration of dilution								13.35	uM	0.20
concentration of stock solution								213.60	uM	4.33
purity-corrected concentration of stock solution								170.88	uM	3.47
Cell solution preparation										
		desired		typical error				actual		error
dilution factor		0.05852		0.00078				0.04676		0.00081
stock solution volume		0.293	mL	0.003	mL			0.234	mL	0.004
buffer volume		4.71	mL	0.05	mL			4.77	mL	0.03
total volume for titrate	(need ca. 2.1 mL/experiment)	5.00	mL	0.05	mL			5.00	mL	0.030
titrate concentration		10.00	uM	0.13	uM			7.99	uM	0.213
Parameters from Origin fit										
		reported from fit		error from fit		percent error from fit		actual		error
n (stoichiometry, purity, and V ₀ correction)		0.628		0.005		0.51%		0.928		0.006
K (association constant)		1.25E+08	M-1	3.99E+04	M-1	3.19%		1.25E+08	M-1	4.24E+04
K _d (dissociation constant)		8.00E-07	M	2.55E-08	M	3.19%		8.00E-07	M	2.71E-08
DeltaH		-10.8	kcal/mol	0.07304	kcal/mol	0.68%		-10.80	kcal/mol	0.14
TDeltaS		-2.48	kcal/mol	0.08	kcal/mol	3.04%		-2.48	kcal/mol	0.14
DeltaG		-8.32	kcal/mol	0.02	kcal/mol	0.23%		-8.32	kcal/mol	0.02